

SYSTEM ENVIRONMENT REQUIREMENTS

Document number: PC-007-System Environment Requirements Rev.15_EN
Document revision date: 2026-01-13

Table of content

1 GENERAL INFORMATION 3

1.1 WARNING3

2 VERIFIED SYSTEM COMBINATIONS 4

2.1 OPERATING SYSTEM4

2.2 NVIDIA GPU4

2.3 NVIDIA GPU DRIVER4

2.4 MINIMUM SYSTEM SPECIFICATIONS5

2.5 REMOTE ACCESS AND VIRTUALISATION PLATFORMS5

2.6 OTHER REQUIREMENTS.....5

3 CONTACT INFORMATION..... 6

3.1 MANUFACTURER CONTACT INFORMATION.....6

3.2 REGULATORY REPRESENTATIVES.....6

3.3 SUBSIDIARIES.....6

1 GENERAL INFORMATION

Hermes Medical Solutions' software has been verified on systems with the specifications described in this document. For questions regarding other combinations, please contact the support service, see Contact Information.

For our reconstruction and voxel-based dosimetry software, there are specific requirements which are described in separate columns and supersedes the general requirements.

Hermes Medical Solutions welcomes feedback from readers of this manual. Please report any errors in content or typography and suggestions for improvements to our support, as noted in Contact Information.

NOTE: It is the customer's responsibility to install and maintain anti-virus software on the server and client computers and apply the necessary protection against potential threats.

1.1 Warning



Only properly trained service personnel of Hermes Medical Solutions, or an authorized partner, shall perform installations, and service of this product.



Using software other than that provided by Hermes Medical Solutions may lead to reduced performance and, in the worst cases, incorrect output data, if installed on the same system as Hermes Medical Solutions software.

2 VERIFIED SYSTEM COMBINATIONS

For all our products - Hermes licensing software: ExplorerPWD V4.6.2 or later.

2.1 Operating System

- Microsoft® Windows 11
- Microsoft® Server 2016
- Microsoft® Server 2019
- Microsoft® Server 2022

2.2 NVIDIA GPU

The installed GPU must be compatible with the corresponding Operating System, GPU Driver, and Hardware/Hypervisor.

NOTE: For GPU-accelerated reconstruction and voxel-based dosimetry at least 6GB GPU Memory is required.

NVIDIA® Series	Validated GPU
Quadro®	Px000/Px200
	Tx000
Quadro® RTX™	x000
RTX™	Axx00
	x000 Ada
Data Center / Tesla	T4
	L2/L20
	L4/L40/L40S

2.3 NVIDIA GPU Driver

The installed GPU Driver must be compatible with the corresponding Operating System, GPU, and Hardware/Hypervisor.

Confirm with NVIDIA that the chosen driver branch is still supported to receive security and bug fixes.

Driver Branch	Minimum Patch
RTX / Quadro Enterprise driver	
R470	U12
R525	U7
R535	U2
R550	U11
R580	U5
Data Center driver	
16	2

2.4 Minimum system specifications

WORKSTATION

- CPU:** Intel® (Recommended 13th Gen Intel Core i7-13700K 16C 3.40GHz or similar)
- RAM:** Recommended minimum 10GB (excluding OS requirement)
- Resolutions:** Tested Min/Max: 1920x1080 up to 3840x2160
Tested: 1920x1080, 1920x1200 and 2560x1440 for Hybrid Recon and Voxel Dosimetry
- Scale:** Use recommended scaling based on resolution
Use 100% for Voxel Dosimetry and Hybrid Recon
- Network:** Recommended minimum 1Gbps

APPLICATION SERVER

- CPU:** Intel® (Recommended Intel® Xeon® Gold 6434 8c 3.7Ghz or similar)
- RAM:** Recommended minimum 10GB per concurrent user (excluding OS requirement)
- Network:** Minimum 1Gbps (10Gbps recommended)

THIN CLIENT (for Application Server access)

- Resolutions:** Tested Min/Max: 1920x1080 up to 3840x2160
Tested: 1920x1080, 1920x1200 and 2560x1440 for Hybrid Recon and Voxel Dosimetry
- Scale:** Use recommended scaling based on resolution
Use 100% for Voxel Dosimetry and Hybrid Reconstruction
- Network:** Minimum 1Gbps (10Gbps recommended)

2.5 Remote access and virtualisation platforms

The software may be deployed in on-premises or remote virtualized environments. The medical device does not rely on or incorporate any functions specific to third-party remote access platforms (e.g. Citrix Workspace, Microsoft RemoteApp, VMware Horizon).

Any remote access or virtualization solution may be used, provided that the healthcare institution ensures that the selected environment meets and maintains:

- acceptable application performance,
- adequate image quality for clinical interpretation,
- compliance with local cybersecurity and IT policies.

Verification and validation of the remote access environment is the responsibility of the healthcare institution.

For guidance on remote access platforms that Hermes Medical Solutions has tested please refer to PC-009 System Environment Guidelines.

2.6 Other requirements

For licensing in cloud environments, Corridor4DM requires a dedicated virtual machine and Cedars-Sinai Cardiac Suite requires either a dedicated virtual machine or outgoing internet access.

3 CONTACT INFORMATION

Contact any of the addresses below for service, support or if you have any other questions.

3.1 Manufacturer contact information

**Head office****Hermes Medical Solutions AB**

Strandbergsgatan 16

112 51 Stockholm

SWEDEN

Tel: +46 (0) 819 03 25

www.hermesmedical.com

General e-mail address:

info@hermesmedical.com

Support e-mail addresses:

support@hermesmedical.com

support.ca@hermesmedical.com

support.us@hermesmedical.com

3.2 Regulatory Representatives

UK Responsible Person

Hermes Medical Solutions Ltd

Cardinal House

46 St. Nicholas Street

Ipswich, IP1 1TT

England, United Kingdom

CH Authorized Representative

CMI-experts

Grellinger Str. 40

4052 Basel

Switzerland

Australian Sponsor

Cyclomedica Australia Pty

Ltd 4/1 The Crescent,

Kingsgrove,

Sydney 2208

Australia

3.3 Subsidiaries

Hermes Medical Solutions Ltd

7-8 Henrietta Street

Covent Garden

London WC2E 8PS, UK

Tel: +44 (0) 20 7839 2513

Hermes Medical Solutions, Inc

2120 E. Fire Tower Rd, #107-197

Greenville, NC 27858

USA

Tel: +1 (866) 437-6372

Hermes Medical Solutions Canada, Inc

1155, René-Lévesque O., Suite 2500

Montréal (QC) H3B 2K4

Canada

Tel: +1 (877) 666-5675

Fax: +1 (514) 288-1430

Hermes Medical Solutions Germany GmbH

Robertstraße 4

48282 Emsdetten

Deutschland

Tel: +46 (0)819 03 25