



Version of SER

H
E
R
M
E
S

P
R
O
D
U
C
T

		Rev.12	Rev.10	Rev.9	Rev.8	Rev.7	Rev.6	Rev.5	Rev.4	Rev.3	Rev.2
				First SER without Server 2012R	Only for MDR reviews	Windows 11, Win 2019 (Citrix / RDP)	Win 2019 (Citrix)	For BJC only	For BJC only	For NA market, Affinity 4.0.0 only	
Affinity	3.0.5					X	X				X
	4.0.0	X	X	X							
	4.0.1	X	X	X							
	4.0.2	X	X					X		X	
	5.0.1	X	X								
Hybrid Viewer	6.1.4					X	X		X		X
	7.0.2	X	X	X	X						
	7.2.0	X	X								
Hybrid Recon	4.0.5								X		
	4.0.6					X	X				X
	5.0.0	X	X	X	X						
Voxel Dosimetry	1.1.0	X	X								
	1.1.2	X	X			X	X		X		
	3.0.0	X	X								
	3.1.0	X	X								
	3.2.0	X	X								
4DM	2018					X					
	2024	X	X	X							
Cedars	2017	X	X	X		X					
Hybrid3D	3.0.2					X	X		X		

Revision number	Revision History	Published
Rev.1	First version of this Matrix. Includes SER versions 2-8. Includes applications: Affinity 3.0.5 and 4.0.0 Hybrid Viewer 6.1.4 and 7.0.0 Hybrid Recon 4.0.5, 4.0.6 and 5.0.0 Voxel Dosimetry 1.1.0, 1.1.2 and 2.0.0 OLINDA/EXM 2.2.3 Hybrid3D 3.2.2 4DM 2018 Cedars 2017 GOLD WS 2.16 and 2.17 GOLD Lx 2.16 and 2.17	1/13/2024
Rev.2	Added tab for EOL and End-of Service dates Added SER Rev.9 Removed Affinity 4.0.0/SER Rev.7, which was incorrect Added 4DM 2024 Changed Voxel Dosimetry 2.0.0 to 3.0.0 Changed Hybrid Viewer 7.0.0 to 7.0.1	3/4/2024
Rev.3	Updated SER Rev.9 to include Affinity 4.0	3/11/2024
Rev.4	Updated SER Rev.9 to include Cedars 2017	3/25/2024
Rev.5	Updated SER Rev.9 to include Hybrid Viewer 7.0.1	4/12/2024
Rev.6	Changed Hybrid Viewer 7.0.1 to Hybrid Viewer 7.0.2	4/17/2024
Rev.7	Updated SER Rev.9 to include Hybrid Recon 5.0.0	4/18/2024
Rev.8	Updated SER Rev.9 to include Affinity 4.0.1	8/30/2024
Rev.9	Added SER Rev.10, removed GOLD WS, GOLD Lx and OLINDA/EXM	9/20/2024

Rev.10	Updated SER Rev.10 to include Affinity 4.0.2	10/30/2024
Rev.11	Updated SER Rev.10 to include Voxel Dosimetry 3.1.0	12/18/2024
Rev.12	Added SER Rev.12 for all products that had SER Rev.10 and for Affinity 5.0.1 Updated SER Rev.10 to include Affinity 5.0.1	5/9/2025
Rev.13	Updated SER Rev.12 to include Hybrid Viewer 7.2.0 Updated SER Rev.10 to include Hybrid Viewer 7.2.0	7/3/2025
Rev.14	Updated SER Rev.12 to include Voxel Dosimetry 3.2.0 Updated SER Rev.10 to include Voxel Dosimetry 3.2.0	9/10/2025